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Establish brand
with quality

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Brief Introduction

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Jiangsu Huatai Heavy Industries Co., Ltd. namely original Danyang Huatai Pharmaceutical Machinery Co., Ltd. is a technology-oriented enterprise which specializes in R & D and manufacture of pharmaceutical machinery and food machinery. With twenty years of experience in production of pharmaceutical machinery, the company owns strong product development, design capacity and rich production experience as well as the capacity of producing third grade pressure vessels by means of complete processing equipments, perfect technological process and strong technical force.

Regarding scientific and technological innovation as the development impetus and core competitiveness, our company has owned many proprietary intellectual property rights and has cooperated with many scientific institutions, universities and research centers, and has absorbed advanced technology and experience from abroad. We make continuous innovation and reform in the research and development of pharmaceutical machinery to have designed and manufactured a series of products with advanced technology. Our company was the first enterprise to produce ultrasonic continual dynamic countercurrent extraction

equipments. In 1992, in order to meet the market demands our company developed continual dynamic countercurrent extraction equipments independently which has won many national patents and reached international advanced level through experts' appraisal. They have been repeatedly recognized as superior pharmaceutical machinery by many domestic companies such as Huabei Pharmaceutical, Changchun Haiwang, Damin Food, Chifeng Danlong, Guilin Sanjin etc. Besides domestic application, these equipments have also been sold far to USA, Uzbekistan and other countries. Our company takes the lead in production of updated products of modernized Chinese herbal medicine, having achieved good social benefits.

Our company adheres to the quality policy of "good faith foremost, scientific management, innovative development and excellent service" and will welcome friends from various industries to visit us for sincere cooperation and common progress by relying on best quality and first-rate service.

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

Honors that have been achieved are due to the dedication of every Huatai member, encourage and urge everyone to seek for high quality and perfect projects. Being unsatisfied with the present and stepping into new stage continuously is the the long lasting faith for Huatai members to forging ahead with determination.

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First-rate quality is embodied by advanced automatic production equipments and superior production line. Direct involvement of technical experts in the production ensures the stability and reliability of products. Steadily increasing awareness of quality products and constant flow of industrial quality products make Huatai Pharmaceutical advance with each passing day. No matter how the world changes, Huatai members will insist on the spirit of seeking for greater perfection continuously.

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Internationalization of equipments is the important guarantee of enhancing international competitiveness.

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We are only dedicated to one industry on a professional basis. Rich experience in domestic and overseas markets as well as deep understanding upon the cultural differences of various countries ensures the unique advantages of Huatai in aspects of design and manufacture, with the ability of offering innovative solutions for you from design to manufacturing process.

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
Process flowchart and small equipment picture

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Contents

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01 超声波逆流提取机组
Ultrasonic Counter-current
Extractor

03 多功能提取罐
Multi-Functional Extractor

05 实验室微型提取浓缩机组
Lab Multi-Functional Extractor and
Evaporator Set

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Purpose

The unit is a kind of brand new dynamic extracting equipment that extracts the effective components in natural products (esp. TCM materials). It has been granted the national patent. Used with other auxiliary equipment such as condenser, it can also meet technological requirements of solvent recovery and etc.

Features

According to modern production theories and principles of natural products (esp. TCM materials), the unit is designed and developed through many years of research and development.

The unit is composed of certain amounts of extraction single tanks provided with special devices and equipped with the heating jacket outside, which connect with material and liquid conveying pipes (the enhanced travel aims to create the maximum strength difference). The extracting temperature can be appropriately selected according to material natures. Generally the steam is adopted as the heat carrier, or hot water or heat transfer oil (to facilitate automatic control under low temperature) can also be adopted. According to TCM extraction theories, generally the optimal temperature during medicine extraction is about 80℃, so as to avoid damaging the heat sensitive components, make starch not easily to be lixiviated and reduce the difficulty of the next filtering procedure.

The unit adopts the continuous countercurrent technology. The medicine materials are fed automatically in fixed quantity from the feeding inlet, the solvent is filled in controlled quantity from the liquid inlet. Under the driving of the special devices, the leaching solution and medicine materials form reverse continuous countercurrent leaching process. Those flowing out from the beginning are the extraction liquids that contain effective components, and those discharged finally are medicine residues that contain no effective components. The residues will enter the extruding machine or solvent recovery machine for disposal and cleaning, which is cost saving and does not pollute the environment. The whole unit is a closed system, which makes it suitable for both extraction of organic solvents and water solvents. It is clean and sanitary. The feeding and residue discharging are automatically completed with convenient operation. During the cleaning, just open the upper door cover. It conforms to GMP requirements.

The equipment has been successfully applied to extraction of TCM, steviol glycosides, tea powder, momordica glycosides, tobacco sheets, etc. With high extracting efficiency and low running cost, it has been received good comment from customers.

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Ultrasonic Counter-current Extractor

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Auxiliary system of ultrasonic continuous countercurrent extracting equipment

1. Medicine materials conveying and control system
2. Ultrasonic generation system
3. Extracting solvent flow and temperature control system
4. Heating and heat preservation system during the extraction process
5. Medicine liquid and residue separation system
6. Medicine residue disposal system: composed of medicine residue extruding machine, drying machine, organic solvent recovery condenser, medicine residue conveyor, etc.
7. Online cleaning system
8. Electric operation and control system: The whole system adopts PLC computer operation and control.

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Purpose

The equipment is applicable to the multiple technology operations of normal pressure, micro pressure, water frying, warm soaking, thermal reflux, compulsory circulation, filtration, redolent oil extraction and organic solvent recycling in traditional Chinese medicine, food and chemical industry. It has the advantages such as high efficiency, simple operation and etc. The CIP cleaning system is set inside the tank, conforming to GMP standard.

Features

1. The residue discharging door adopts ordinary double-cylinder starting and closing type or three-cylinder rotary safe door type. The double-cylinder rotate to push the clamp ring for fastening and the bevel wedge locks itself, so as to solve leakage or disengagement caused by unstable pressure of the compressed air supply and realize safe production.
2. It can be applied in the dynamic extraction system after the liquid feeding mode is changed and the blending device is added.
3. The annular filter screen is additionally provided at the bottom of the tank to enhance the area of the filter screen. The area of the filter screen is several times higher than that of the old style extraction tank.

4. The pump jetting forms upward forced countercurrent circulation, playing the function of blending.
5. The new type foam collector is adopted to prevent the blocking formed by the small and light medicine materials in the foam that enter the cooler lined tube.

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Lab Multi-Functional Extractor and Evaporator Set

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
Purpose:

The equipment is suitable for extraction of TCM and various natural plants. It is the ideal device for determination of technological parameters before modernized industrial production

The whole set of equipment has complete functions, and the whole system vacuum can be achieved. It can meet various extraction requirements such as water extraction, alcohol extraction, extraction in batches, warm soaking and heat reflux extraction, percolation extraction, etc. Extraction, chromatography, filtration, concentration, drying, ointment collection, organic solvent recovery, etc. can be carried out in the extraction solution. The optimal combination of various technologies can be realized. The system has two modes such as steam heating and electric heating for users to select, so it is especially suitable for units without steam.

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Purpose:

The Soxhlet extracting tank conforms to GMP standard and has good technological adaptability. It is suitable for water extraction and alcohol extraction under conditions of normal pressure, negative pressure and positive pressure, especially suitable for low-temperature extraction and concentration of heat sensitive materials. Low-temperature water extraction can be carried out with temperature above 45℃. The specific gravity of concentration liquid can reach 1.1-1.35. For special materials (likely to agglomerate and solidify), natural cycling can be changed into compelled cycling to meet customer requirements.

Highly automatic degree. The temperature, pressure, flow, liquid level and concentration can be set and kept under automatic control with convenient operation and stable performance. All instruments, meters, executive elements and PLC are made of famous brand products from home and abroad. Therefore, they are highly reliable.

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Purpose:

This equipment is widely used in medicine, light industry, food industry and the evaporation of liquid material mixing process needs income, especially for low-temperature vacuum heat-sensitive materials, vacuum concentration, can produce up to 1.4 weight extract. Patent number: ZL92218691.X.

Features:

- 1, using two effect evaporator while making full use of secondary steam. evaporation chamber to reduce the evaporation temperature for heat-sensitive drugs
- 2, 2-effect evaporation chamber wall by scraping device, to prevent the process of coking, enrichment up to 1.4 specific gravity
- 3, PLC control, and achieve automatic condensate discharge
- 4, by external heating cycle of intense heating process, evaporation speed, is short
- 5, all of high quality stainless steel equipment, comply with food hygiene of Traditional Chinese Medicine, to GMP requirements

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Purpose:

This concentrator is applicable to the concentration for liquid material of Chinese traditional medicine, Western medicine, glurmet powder, dairy food, chemical engineering, especially suitable for concentrate of heat sensitive materials in vacuum under low temperature. It has more advantages than previous products after careful calculation and elaborative improvements in design. Appropriate concentrator can be selected in accordance with its concentration capacity.

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Purpose:

This equipment is applicable to the concentration for materials of Chinese traditional medicine, starch sugar etc. and recovery for industrial organic solvent (like alcohol), and is suitable for low temperature vacuum concentration of heat sensitive materials with small quantity and great category.

Alcohol recovery: recovery capacity, vacuum concentration process. Than the old type of similar equipment to improve productivity 5-10 times and reduce energy consumption 30%, with investment, recovery efficiency is high.

Concentrate feed: This equipment is used outside the natural cycle of heating and vacuum evaporation of a combination, evaporation speed, concentration ratio up to 1.3.

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Ultralow-temperature Rotary Film Evaporator

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Purpose:


The equipment is suitable for plant extract, traditional Chinese medicine, biological products, fruit juice, milk products, instant tea, which needs to be concentrated and recovered under Ultra-low temperature conditions, and requires that the equipment can be cleaned thoroughly and conveniently

Features:

1. Continuous feeding and continuous discharging under vacuum condition;
2. Evaporate process achieved on the pipeline, continuous, automatic closed operation;
3. Ultra-low temperature (30 -50 C) can be adjusted to ensure that the heat sensitive ingredients are not damaged;
4. According to the characteristics of the material, the evaporator can be CIP online cleaning, and can also move the heater outside for cleaning, no dead corners, blind spots, in line with the requirements of GMP specification;
5. Recoverable solvent in operation;
6. It also has the low altitude model adapted to the low workshop

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Product Usage:
Pure electric concentrator imitates the refrigeration principle and uses a carrier medium to repeatedly circulate gas and liquid in a closed cycle system,thereby achieving heat exchange for heat release and absorption. During this process,the released heat energy is used to heat the liquid and promote its evaporation,while the absorbed heat energy comes from the condensation of secondary steam generated during the evaporation of the liquid.Through the dual utilization of cold and hot energy,pure electric concentrator can achieve the concentration of feed liquid,evaporation of solvent,and recovery of solvent.Compared with traditional heating methods, the energy consumption is significantly reduced,the temperature control accuracy is high,the degree of automation is high,the modular design is easy to maintain,and the safety performance is high.

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Purpose:

OZR series centrifugal granulator dryer suitable for ceramic industry, pharmaceutical industry, chemical industry, dye and other special materials, centrifugal granulation, particle size distribution requirements can be adjusted to the optimum distribution curve of values (100 mesh -350 mesh), its The main technical performance, indicators can reach the world advanced level.

Features:

Fast drying, material after atomization surface area increased significantly in the hot air flow, the moment could evaporate 95% -98% of the water, only a few seconds to complete the drying time, especially suitable for drying heat sensitive materials. Products with good uniformity, mobility, and solubility, high purity and good quality. Production process is simplified, convenient control. 40-60% moisture content for special materials, 90% of liquid can be dried into a powder product, after drying without crushing and screening to reduce the production processes, improve product purity. Product particle size, bulk density, moisture, in a certain range can be adjusted by changing the operating conditions, control and management are very convenient.

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Purpose:

The equipment used for pharmaceutical and chemical industry, the recovery of the purification of dilute alcohol, also apply to methanol, acetone and other organic solvent distillation recovery. Features: A device operation method of the two kinds of equipment can be 30 °C - 50 °C dilute alcohol distillation to 90 °C - 95 °C residue emission-grade = 3 °C, is the preferred alcohol recovery specialist equipment.

Features:

This device adopts high efficiency stainless steel corrugated padding to ensure finished alcohol will not be polluted. Tower still adopts detachable U-shaped tube, which equipment can adopt two modes, external feedback heating and internal feedback heating and users can purchase what they need according to actual demands.

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Purpose:

Hot air circle oven is used in the pharmaceutical, chemical, food, agricultural and sideline products, aquatic products, light industry, heavy industry materials and products such as heat curing, drying dehydration. Such as raw materials, crude drugs, herbs and herbal extract, powder, granule, powder, water pills, bottles, pigment dyes, dehydrated vegetables, dried fruit, sausages, plastic resin, electric elements, such as baking finish. This oven features the largest series in the box most of the hot air to circulate, thus increasing the difference between the upper and lower wet. The equipment available in computer-controlled system to realize automatic temperature control and fully enclosed, full-cycle operation.

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Purpose:

This machine is suitable for pharmaceutical, chemical, food, electronics, medicine and other industries drying.

Features:

- 1.Vacuum lower boiling point material solution to the evaporator heat transfer driving force increases, and therefore a certain quantity of heat, can save the evaporator heat transfer area.
- 2.Applicable to high temperatures that can decompose and polymerize and metamorphose thermal drying process
- 3.During the drying process were not mixed with any impurities, material and physical harm.

技术参数 Technical parameters

项目\型号 Item \ Model	YZG-600	YZG-1000	YZG-1400	FZG-15
干燥箱内尺寸mm Drying in size	976	1527	2054	1500x1400x1220
干燥箱外尺寸m Oven outside dimensions	1.14x0.81x1.02	1.7x1.2x1.5	2.4x1.67x1.9	1.5x1.9x2.11
烘架层数 Layer drying rack	4	6	8	8
层间距离mm Interlayer distance	81	102	102	102
烘盘尺寸mm Baking tray size	310x600x45	250x410x45	400x600x45	480x630x45
烘盘数 Number of baking trays	4	24	32	32
干燥箱重量kg Drying weight	250	800	1400	2100

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Vacuum Decompression Concentrator

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Purpose:

This equipment is mainly used in pharmaceutical, food, health care products, oral liquid, juice, chemical and other industrial materials, concentrated liquid evaporation, and also recovery of alcohol and simple back extraction.

Features:

This equipment includes concentrated solution, the first condenser, gas-liquid separator, cooler, bucket six parts by the fluid, all stainless steel tank for the concentrated jacket structure for the column tube condenser, cooler for the coil type.

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Purpose:

Is a heat exchanger to achieve the transfer of materials between the heat energy-saving equipment, in the petroleum, chemical, petrochemical, metallurgical, power, light industry, food industry, such as the universal application of a technology equipment. In oil refineries, chemical installations in the heat exchanger of the total number of equipment about 40 percent of the total investment of 30% -45%. In recent years, with the development of energy-saving technologies, applications continue to expand, using a high-temperature heat exchanger and low-temperature heat recovery has brought significant economic benefits

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
Purpose:

This machine is suitable for all kinds of herbs, woody plants, residues of residue to squeeze out the juice, extract recovery solvent and the active ingredients to enhance the leaching rate of more products, increase cost effectiveness. This machine is applied to all pharmaceutical companies, research institutes, health care products with herbs and vegetable plants and fruit as raw materials for food, beverage, cosmetics, make wine and other industries.

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Purpose:

This equipment can be used for medicine, food, chemical industry, the liquid material storage. High efficiency, convenient operation, the tank is equipped with auto-rotating spray head to ensure a complete cleaning. Materials are stainless steel SUS304 or SUS316L manufacturing, internal surface polishing to $Ra \leq 0.4\mu m$, external surface matt treated to meet the GMP pharmaceutical standards. Stainless steel plate system for tank cars, horizontal, vertical tank, high chamfer. Its, special size according to user requirements design and manufacture.

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Purpose:

This equipment is applicable to alcohol precipitation in Chinese medicine, oral solution, food, health products, chemical engineering etc. Alcohol precipitation tank is kind of newly-developed special-purpose equipment by adopting old technology, namely, depositing. It is mainly applied to cooling operation of Chinese medicine decoction after concentration or alcohol precipitation, and can also be used for water precipitation operation of concentrated solution of extracted Chinese medicine.

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Purpose:

This mixing tank is also called batching bowl, material stirring bowl, batching barrel, mixing tank, with advantages of energy-saving corrosion resistance, strong production capacity, easy cleaning, simple structure etc. It is mainly used to evenly mix dairy product and sugar as well as other elements and carious materials after they are mixed. It is a kind of indispensable equipment for manufacturers of dairy products, beverage and medicine.

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Purpose:

The equipment used in pharmaceutical, chemical, food and light industry in the hydrolysis, neutralization, crystallization, distillation, evaporation and other processes.

Features:

Mixing patterns have anchored, frame type, blade type, etc.; rotating body can use cycloidal reducer, stepless speed, gear, etc.; seal mechanical seal can be used; heating and cooling jacket can be used, half pipes, plates tube structure; heating methods are steam, electric heating, oil, etc.; to meet the acid and alkali, high temperature, wear, corrosion and other processes in different working environment environment technology needs. Process according to user requirements can also design and manufacture.

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Purpose:

Stainless steel pan, also known as sandwich sandwich steam pot, widely used in confectionery, pharmaceutical, dairy, wine, cakes, preserves, beverages, canned, marinated and other food processing, food processing to improve quality, reduce time, improve working conditions good equipment.

Features:

High efficiency, uniform heating, easy control, can effectively shorten the boiling time of the material liquid, beautiful and practical.

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Purpose:

This machine is suitable for chemical, pharmaceutical, food and other industries powder, satellites and fiber enrichment materials, mixing, drying and the need for low-temperature drying of the material (such as a high sugar content, viscosity of the thick material), is more suitable the easily oxidized, volatile, heat-sensitive materials are dry and dried granule of proprietary Chinese medicines. The machine used to seal mechanical seal, can avoid the run, run, drain, drop, and frequently changed filler Zhi Wei Xiu hard to conform to GMP on equipment Yaoqiu conditions, the machine has insulation layer, Nei Wai material system 1CR18Ni9Ti stainless steel, simple, easy to clean.

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Jiangsu Huatai Heavy Industries Co.,Ltd.